

Document No. 3104, A03, A02	Revision Number 00	Effective Date 07 July 2008	Supersedes Date N/A	Page 1 of 1
Document Title: PURIFIED HUMAN PANCREATIC ISLETS, GLUCOSE STIMULATED INSULIN RELEASE DETERMINATION BY ELISA, ATTACHMENT II, STATIC INCUBATION & INSULIN RELEASE WORKSHEET				

Follow instructions in SOP 3104, A03

	Pre-culture		Product Bag #1		Product Bag #2		Product Bag #3	
	Low Glucose	High Glucose	Low Glucose	High Glucose	Low Glucose	High Glucose	Low Glucose	High Glucose
Stimulation								
Start Time								
Stop Time								
Stimulation Time								
Performed by								
Date								

Enzyme-conjugate Incubation	Pre-culture		Product Bag #1		Product Bag #2		Product Bag #3	
Start Time								
Stop Time								
Incubation Time								
Shaker RPM								
Room Temperature								
Substrate TMB Incubation								
Start Time								
Stop Time								
Incubation Time								
	Low Glucose	High Glucose	Low Glucose	High Glucose	Low Glucose	High Glucose	Low Glucose	High Glucose
Optical Density								
mU/L								
Stimulation Index								
Performed by								
Date								

Reviewed by: _____ Date: _____

Pre-culture Islets Lot #: _____ Islets Product Lot #: _____